



**Imaging Spinal Surgical Tabletop
Models: 5827, 5927
of the Modular Table System (MTS) product line
Work Instruction Number 50603-12**

Contact Information:

Date: _____	Model/Serial#: _____	Account #: _____
Hospital Name: _____		City & State: _____
Contact Name: _____	Position: _____	Phone: _____
Equipment Location: _____		
Service Representative: _____		Service Organization: _____
Service Rep. phone: _____		Home location: _____
Customer contact accepting service: _____		
Name please print		Signature
		Date
Service Representative: _____		
Name please print		Signature
		Date

Using this report: Please proceed through each step and place an "x" in the "()" as each item representing Pass Fail () () criteria is inspected. If a step has failed, please indicate with an "x" and note the reasons in the comments section. If you resolve the failure place an "x" in the "Pass" area and indicate the corrective action again in the comments section. If you are unable to perform any Pass/Fail check, please contact Mizuhosi Technical Services for assistance and circle the check.

Biohazards when servicing equipment: Please take appropriate precautions for servicing potentially contaminated biohazard equipment or working in biohazard surroundings. Please wear protective surgical gloves and log sleeved gowns at all times. If your gloves or protective clothing tears, replace immediately. Have the hospital staff clean the equipment appropriately prior to servicing. Clean your tools that may contain biohazards prior to placing back in your tool kit for removal from the location. If you suspect biohazards are present do not service the equipment until they have been cleaned.

(If serviced by an outside service organization, the record must be reviewed, signed, and dated by Technical Services Dept. Refer to last page of this document).

Surgical bases covered:

- 5827 Radiolucent Imaging Top, 350lb (pounds) patient weight capacity
- 5927 Radiolucent Imaging Top, 500lb patient weight capacity

Tools Required:

- Cotton swab
- Lubricant, silicon, spray
- Grease, Lithium bases, white, multi-purpose

Materials Required:

- Jackson Spinal Surgery & Imaging Table Tops Technique Manual

Pass Fail

1. Top Inspection and fit.

- 1.1 ☐ ☐ Table top is clean, without signs of hazardous contamination, and all labels are intact and legible.

Comment: _____

- 1.2 ☐ ☐ Carbon fiber top is free of wear, cracks, delamination or splinters. Particular inspection is required of the siderail-to-top joint. Joint is securely bonded without voids, cracks, or separation.

Comments: _____

Diagnostic/repair tip: splinters are easily detected by running a cotton swab along surface. Smooth splinter area with a *light* sanding with 400 grit sandpaper then applying a *small* amount of Polycyanuacrylate ("super-glue") to seal affected area.

- 1.3 ☐ ☐ Gimbals are secure, rotate to stops and return to center position.

Comments: _____

Note: earlier models may not be equipped with gimbals.

- 1.4 ☐ ☐ Fits in all positions on removable H-frames, T-pins slide into locked position freely.

Comments: _____

Warning: Rotation of tops not equipped with gimbals beyond 5 degree Trendelenburg and Reverse-Trendelenburg will cause severe damage.

- 1.5 ☐ ☐ With top in 3rd from highest position on H-frames, top rotates freely at 5 degree Trendelenburg and Reverse-Trendelenburg.

Comments: _____

- 1.6 ☐ ☐ With top in lowest position on H-frames, top rotates freely at 5 degree Trendelenburg and Reverse-Trendelenburg.

Comments: _____

2. Standard Attachments operate properly.

2.1 ☐ ☐ Table pad retains good shape, is free of punctures, tears or excessive wear.

Comments: _____

3. Optional Accessories inspection only if they are attached to the top. Otherwise, no need to inspect.

3.1 ☐ ☐ Pivoting Armboard(s) operates properly and the pads is/are free of punctures, tears or open seams.

Comments: _____

3.2 ☐ ☐ Clarke Socket(s) operate properly.

Comments: _____

3.3 ☐ ☐ Cross Arm Support(s) operate properly and the pads is/are free of punctures, tears or open seams.

Comments: _____

3.4 ☐ ☐ Aluminum Siderail Socket(s) operate properly.

Comments: _____

3.5 ☐ ☐ Imaging Top Adaptor Bracket(s) operate properly.

Comments: _____

3.6 ☐ ☐ Siderail Adaptor(s) operate properly.

Comments: _____

This section is for Technical Services use only.

Reviewed: _____

Signature: _____

Date: _____